



ROTEX™ 780

COMPACT ELECTRIC HANDPIECE UNIT

ROTEX™ 780

- Ball bearing high torque electric motor.
- 3,500 to 35,000 R.P.M. variable speed.
- Double insulation heavy duty construction.
- Automatic cooling fan for continuous high speed performance.
- Forward/Reverse switch.
- Standard package includes Foot Switch, Handpiece, and Handpiece Rest.
- Optional variable speed foot control.
- Application for grinding, polishing, milling of precious metals, ceramics, acrylics, and steel.
- Unit is available for 100/110V 50/60Hz, or 220/240V 50Hz operation.

ISO 9001
ISO 13485



DENTAMERICA®

18320 Bedford Circle
P.O. Box 3200
Industry, CA 91744 U.S.A.

Phone: (626) 912-1388
Fax: (626) 913-0510
www.dentamerica.com